



Diacclone

DIACLONE Immunology Products

IL-10 ELISA Set

Human IL-10 ELISA Set

CATALOGUE N° :	851.540.001, 851.540.005, 851.540.010, 851.540.015, 851.540.020
SPECIFICITY :	Recognizes both natural and recombinant human IL-10
RANGE :	Recommended from 12.5 pg / ml to 400 pg / ml
SENSITIVITY :	< 5 pg / ml
INCUBATION :	From sample to end 2h45
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human cytokines Cross reactivity with simian IL-10 No cross reactivity with viral IL-10
KIT SIZE :	1X96 (Discovery), 5X96, 10X96, 15X96, 20X96
REFERENCES :	Audran, R. et al., Infect Immun., 2005; 73(12): 8017-26. Audran, R. et al., PLoS One, 2009; 4(10): e7304 Boniface, K. et al., Clin Exp Immunol., 2007; 150(3): 407-15. Carrion, M. et al., J. Leukoc. Biol., 2016; 100(6): 1385-1393 Filaci, G. et al., J. Immunol., 2007; 179(7): 4323-4334 Ghannam, S. et al., J. Immunol., 2010; 185(1): 302-312 Haller D. et al., Clin. Diagn. Lab. Immunol., 2002; 9(3): 649 - 657 Heemann, C. et al., Clin. Cancer Res., 2012; 18 (13): 3637-3647 Popov, A. et al., J. Immunol., 2008; 181(7): 4976-4988. Rask, C. et al., Clin Exp Immunol., 2013 172(2): 321-32 Szpakowski, P. et al., J Immunol Res., 2015: 359153 Teitelbaum, D. et al., Mult Scler., 2003; 9(6): 592-9. Urzainqui, A. et al., J. Immunol., 2007; 179(11): 7457-7465. van Hooij, A. et al., Sci Rep., 2016; 6: 34260. Zhang, Q. et al., Infect Immun., 2006; 74(8): 4735-43.
KIT CONTENT :	Capture antibody, Biotinylated detection antibody, Standard, Streptavidin-HRP, TMB

Version 10 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
– BP 1985
25020 Besançon cedex – France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com