



Diacclone

DIACLONE Immunology Products

TNF- α ELISpot Set

Human TNF- α ELISpot - Enzymatic

| | |
|-------------------------|--|
| CATALOGUE N° : | 856.041.001, 856.041.005, 856.041.010, 856.041.015, 856.041.020 |
| SPECIFICITY : | Recognizes natural human TNF- α |
| INCUBATION : | 3 h after cell stimulation |
| CROSS REACTION : | No cross reactivity with other human cytokines Cross reactivity with macac TNF- α No interference with TNFR I and TNFR II |
| KIT SIZE : | 1x96 (Discovery), 5X96, 10X96, 15X96, 20X96 tests |
| REFERENCES : | Godard B. et al., Hum. Immunol., 2004; 65(11): 1307-18 Hudak, S. et al., J Immunol., 2002; 169(3): 1189-96. Kalogerakou, F. et al., Hippokratia, 2008; 12(4): 230-5. Capture and detection antibodies, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer. Kit can be sold with or without plates. |
| KIT CONTENT : | To receive kit plus non-sterile plates (MSIPN), please add 'P' at the end of the Cat N°. To receive kit plus sterile plates (MSIPS), please add 'S' at the end of the Cat N°. |

Spots read visually or with reader

Version 8 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985

25020 Besançon cedex - France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

info@diacclone.com



www.diacclone.com