

MycoGenie Mycoplasma Detection Kit

Visual Mycoplasma Detection



Visual Mycoplasma Detection

Detect Mycoplasma with a simple visual color change!

Rapid Cell Culture Results

1-step, 60-minute assay with no special instrumentation required

Highly Sensitive Assay

Detect as little as 500cfu Mycoplasma in 1ul of cell culture supernatant


Widely Validated

Detect up to 28 Mycoplasma species that contaminate cell cultures

Cost-Effective

Low-cost alternative to PCR, QPCR & ELISA detection methods

Visual detection of Mycoplasma contamination



Tube 1: Positive Mycoplasma Control

Tube 2: Day 0

Tube 3: Mycoplasma Eliminated Sample

Tube 4: Negative Control

Discount


CODE	PRODUCT	PACK SIZE	PRICE	PROMO PRICE
MORVOO11	MycoGenie Rapid MycoPlasma Detection Kit	20 Tests	299 €	209,30 €
MORVO011-50	MycoGenie Rapid MycoPlasma Detection Kit	50 Tests	427 €	298,90 €












Organoid Cell Viability



YOU HAVE THE VISION.

WE HAVE THE SUBSTANCE.

Apical-Out Airway Organoids were treated with A) 0.1% DMSO Control or with B) 10% DMSO to induce cell death. Cell viability was quantified using the 3D Organoid Viability Assay or a Competitor Product after a 1-hour incubation.

CODE	PRODUCT	PACK SIZE	PRICE	PROMO PRICE
AKES081-100	3D Organoid Cell Viability Assay (100 Assays)	100 Assays	203 €	142,10 €
AKES081-500	3D Organoid Cell Viability Assay (500 Assays)	500 Assays	417 €	291,90 €

